

# Expanding the Practice of Peripheral Nerve Stimulation with a Novel, Battery-free Micro-implantable Neurostimulation System

**Sean Li, MD**

Premier Pain Centers,  
Affiliate of National Spine and Pain Centers  
Shrewsbury, NJ

Adjunct Clinical Associate Professor,  
Rutgers New Jersey Medical School, Newark, NJ  
President, New Jersey Society of Interventional Pain Physicians

Smartphone based \  
Remote Control

Therapy Disc  
(External Power)

Leads

Battery-free, micro-Implantable  
Pulse Generator (mIPG)



**AMERICAN SOCIETY OF  
INTERVENTIONAL PAIN PHYSICIANS**  
THE VOICE OF INTERVENTIONAL PAIN MANAGEMENT

# Disclosure

- Consultant/Independent Contractor: Abbott, Avanos, Averitas Pharm, Biotronik, Boston Scientific, Nalu, Nevro, PainTeq, Saluda, SPR Therapeutics, Vertos Medical
- Grant/Research Support: Avanos, Biotronik, Boston Scientific, SGX Medical, Nalu, PainTeq, Saluda, SPR Therapeutics
- Advisory Board: Abbott, Saluda
- Equity: Nalu, National Spine and Pain Centers


# Underpenetrated patient populations

## Estimated Patient Population

**~13–16m**



estimated US adults in 2019 with chronic, moderate to severe, neuropathic pain

**Potential SCS Patient Population**  
Typically for **chronic neuropathic low back pain (CLBP)**



Estimated adults in the  
United States suffering from  
moderate to severe pain

**Potential PNS Patient Population**  
Typically for **pain of a peripheral origin**



Estimated patients seeing  
specialists





- Underserved by existing products

1. Nalu Medical S-1 SEC Filing: <https://www.sec.gov/Archives/edgar/data/1647424/000119312521211837/d102511ds1.htm>


# Anatomical sites for chronic pain of a peripheral origin

## Potential PNS Patient Population

Typically for pain of a peripheral origin



# Limitations w/Certain Existing PNS Technologies




**Limited  
capabilities**

**Limited system  
configurations**




**No trial**

**Limited  
longevity  
(<60 days)**

# Microelectronics & Advanced Stimulation



# micro-IPG – Bringing SCS capabilities to PNS



## SCS systems

- IPG
- Broad capabilities
- Multiple lead configurations
- Full trial
- Impedance measurements
- Longevity

## PNS systems


- Small size
- Can be placed near target nerve







# micro-IPG - Peripheral Nerve Stimulation Examples

## Internal micro-IPG




## External wearable



These images are provided for illustrative purposes only.  
Exact placement of the IPG and Therapy Disc are determined independently by  
physicians.



# Portfolio of Leads & Micro-IPGs



# Menu of therapy options offered in the micro-IPG

**Conventional therapies**  
(available today)

Tonic

Pulse dosing (high and low dose)

Current steering (SCS)

**Advanced therapies**  
(available today)

Burst (multiple variations)

Scheduling

Paired Therapies

Combination Therapies

**micro-IPG therapies**

Current Steering (PNS)

Ultra-long Pulse Width (up to 2000  $\mu$ s)

Pulsed stimulation patterns (PSP)


**Potential for Software updates**  
subject to FDA review and 510(k) clearance



# Expanding Options with PNS

# The PNS patient – patient identification

- **Patient preference – small incision, no implantable battery**
  - Reluctance to undergo “surgery”
  - Fear of having implanted device and battery
  - Pocket pain
  - Need for battery replacement
  - Cosmetic outcome
- **Clinical factors**
  - History of spine surgeries
  - Co-morbidities
  - Anti-coagulants
- **Technical Implementation of a micro-IPG**
  - Dual lead capable, multi-contact leads w/current steering
  - Trial capability (up to 30 days)
  - Easily upgradeable stimulation capabilities (SCS capabilities in a PNS form factor)



# Ultrasound Targets

## **Lower Extremity**

- Femoral
- Sciatic
- Saphenous
- Tibial
- Common Peroneal

## **Truncal Area**

- Ilioinguinal
- Iliohypogastric
- LFCN


## **Upper Extremity**

- Ulnar
- Radial
- Medial
- Axillary
- Suprascapular




# Fluoroscopy Targets

- Suprascapular
- Axillary
- Superior Genicular
- Intercostal
- Cluneal




# 84F, Chronic knee pain

PNS: left superior medial & lateral genicular nerves






# PNS for Failed Back Surgery Syndrome




- 68-year-old male, “Joseph”
- Chronic low back pain
- Prior laminectomy
- Dx: failed back surgery syndrome
- Failed traditional SCS trial
- 80% pain relief with 60-day PNS
- Indirect micro-IPG PNS implantation (#2 in U.S.)



# PNS Patients



# Not All SCS Candidate Undergo Implant



“The most common reasons for rejecting the trial implantation were either **fear of the stimulator** or **unwillingness to undergo another surgery.**”





- PNS candidates exist today within your practices.
- PNS Today:
  - improved technologies
  - cost effective
  - low risk
  - favorable reimbursement
  - Optimal salvage option (failed 60-day, SCS, DRG, etc.)
- Low barrier of entry - certain PNS targets can be accessed under fluoroscopic guidance